



*Enclosed are Origami Protective Face Shields constructed of 5mil polyester and a non-latex elastic(headband). They were created in a non-sterile, sanitary environment by volunteers at PDX Hackerspace (<https://pdxhackerspace.org>) in Portland, Oregon. You can find more information about these devices here: <https://kik.to/wF>, Disinfection Guidance: <https://kik.to/aX> Assembly Video: <https://kik.to/Ta> Text Instructions: <https://kik.to/Xw>*

## **Prerequisites**

These face shields require wiping (to remove any films from the manufacturing process) and sterilization before usage.

## **Usage**

These face shields are intended for use by medical professionals under this FDA policy advisory: <https://www.fda.gov/media/136449/download>

## **Disclaimer**

THESE DEVICES ARE PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR CREATORS OR RELATED VOLUNTEERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THESE DEVICES OR THEIR USE OR OTHER DEALINGS.

Digital Shipping Note: